

FC5816/FC5816R		0.2/0.4 ml	0.5 ml	1.5/2.0 ml	4.5 ml	5 ml	7 ml	10 ml	15 ml	16 ml	25 ml	30 ml	50 ml	85 ml	100 ml	175/225 ml	250 ml	4-10 ml*	5-7 ml*	PCR	MTP	
Swing out rotor, 4 x 250ml, ID																						
Item	30314820																					
For	FC5816			▶					▶	▶					▶							
Max Speed	4500 rpm																					
Max RCF	3780 x g																					
	3780 x g																					
	3780 x g																					
Swing out rotor, 2 x 3MTP, ID																						
Item	30314824																					
For	FC5816			▶																		
Max Speed	4500 rpm																					
Max RCF	2716 x g																					
	2716 x g																					
	2716 x g																					
Swing out rotor, 16 x 50 ml, ID																						
Item	30314828																					
For	FC5816								▶				●					▶	▶			
Max Speed	4500 rpm																					
Max RCF	3735 x g																					
	3735 x g																					
	3735 x g																					
Angle rotor, 6 x 250 ml, ID																						
Item	30314821																					
For	FC5816			▶					▶	▶												
Max Speed	8000 rpm																					
Max RCF	10016 x g																					
	10016 x g																					
	10016 x g																					
Angle rotor, 6 x 85 ml, ID																						
Item	30314826																					
For	FC5816								▶	▶				▶								
Max Speed	11000 rpm																					
Max RCF	15555 x g																					
	21726 x g																					
	21726 x g																					
Angle rotor, 4 x 85 ml, ID, Hi																						
Item	30314827																					
For	FC5816								▶					▶								
Max Speed	12000 rpm																					
Max RCF	14809 x g																					
	14809 x g																					
	14809 x g																					
Angle rotor, 10 x 50 ml, ID																						
Item	30314829																					
For	FC5816								▶					▶								
Max Speed	9000 rpm																					
Max RCF	11771 x g																					
	16022 x g																					
	16022 x g																					
Angle rotor, 6 x 50 ml, ID, Hi																						
Item	30314831																					
For	FC5816								▶													
Max Speed	13000 rpm																					
Max RCF	15869 x g																					
	15869 x g																					
	15869 x g																					
Angle rotor, 30 x 15 ml, ID																						
Item	30314832																					
For	FC5718			▶																		
Max Speed	4500 rpm																					
Max RCF	2830 x g																					
	2830 x g																					
	2830 x g																					

* Blood collection tube

ID: auto rotor identification

HI: version with higher speed

● Adapter not required

▶ With adapter

(Continue on next page)

Distributed by:
 Safety Emporium
 PO Box 1D03
 Blackwood, NJ 08012
 Ph: (866) 326-5412 toll-free
 Fax: (856) 553-6154
 esupport@safetyemporium.com
 www.safetyemporium.com

FC5816/FC5816R		0.2/ 0.4 ml	0.5 ml	1.5/ 2.0 ml	4.5 ml	5 ml	7 ml	10 ml	15 ml	16 ml	25 ml	30 ml	50 ml	85 ml	100 ml	175/ 225 ml	250 ml	4- 10 ml*	5- 7 ml*	PCR	MTP
Angle rotor, 20 x 10 ml, ID, Hi																					
Item	30314833																				
For	FC5816																				
Max Speed	12000 rpm																				
Max RCF	15775 x g																				
Angle rotor, 44 x 1.5/2.0 ml, ID																					
Item	30314835 (30372718)																				
For	FC5816	▶	▶	●																	
Max Speed	15000 rpm																				
Max RCF	21128 x g																				
Angle rotor, 24 x 1.5/2.0 ml, ID																					
Item	30314838																				
For	FC5816	▶	▶	●																	
Max Speed	15000 rpm																				
Max RCF	21379 x g																				
Angle rotor, 24 x 1.5/2.0 ml, ID, BIOSEALS																					
Item	30304361																				
For	FC5816	▶	▶	●																	
Max Speed	15000 rpm																				
Max RCF	21379 x g																				

* Blood collection tube
 ID: auto rotor identification
 Hi: version with higher speed

● Adapter not required ▶ With adapter