

FC5718/FC5718R			0.2/ 0.4 ml	0.5 ml	1.5/ 2.0 ml	4.5 ml	5 ml	7 ml	10 ml	15 ml	16 ml	25 ml	30 ml	50 ml	85 ml	100 ml	4- 10 ml*	5- 7 ml*	PCR	MTP
<b>Swing out rotor, 4 × 100 ml, ID</b>																				
Item	30314822																			
For	FC5718	FC5718R			▲	▲				▲		▲	▲	▲	▲	●	▲	▲		
Max Speed	4500 rpm	4500 rpm																		
Max RCF	3350 × g	3350 × g																		
<b>Swing out rotor, 4 × 100 ml, ID, Sealable</b>																				
Item	30314823																			
For	FC5718	FC5718R			▲				▲	▲		▲	▲	▲	▲	●		▲		
Max Speed	5000 rpm	5000 rpm																		
Max RCF	3885 × g	3885 × g																		
<b>Swing out rotor, 2 × 3MTP, ID</b>																				
Item	30314824																			
For	FC5718	FC5718R			▲															●
Max Speed	4500 rpm	4500 rpm																		
Max RCF	2716 × g	2716 × g																		
<b>Angle rotor, 6 × 85 ml, ID, Hi</b>																				
Item	30314825																			
For	FC5718	FC5718R								▲	▲		▲	▲	●					
Max Speed	11000 rpm	13500 rpm																		
Max RCF	13932 × g	20984 × g																		
<b>Angle rotor, 6 × 85 ml, ID</b>																				
Item	30314826																			
For	FC5718	FC5718R								▲	▲		▲	▲	●					
Max Speed	9000 rpm	9000 rpm																		
Max RCF	10413 × g	10413 × g																		
<b>Angle rotor, 4 × 85 ml, ID, Hi</b>																				
Item	30314827																			
For	FC5718	FC5718R								▲	▲		▲	▲	●					
Max Speed	12000 rpm	12000 rpm																		
Max RCF	14809 × g	14809 × g																		
<b>Angle rotor, 10 × 50 ml, ID</b>																				
Item	30314829																			
For	FC5718	FC5718R			▲		▲	▲		▲					●					
Max Speed	7500 rpm	7500 rpm																		
Max RCF	8174 × g	8174 × g																		
<b>Angle rotor, 6 × 50 ml, ID</b>																				
Item	30314830																			
For	FC5718	FC5718R			▲		▲	▲		▲	▲		▲		●					
Max Speed	6000 rpm	6000 rpm																		
Max RCF	4427 × g	4427 × g																		
<b>Angle rotor, 6 × 50 ml, ID, Hi</b>																				
Item	30314831																			
For	FC5718	FC5718R								▲	▲		▲		●					
Max Speed	12000 rpm	12000 rpm																		
Max RCF	13522 × g	13522 × g																		
<b>Angle rotor, 30 × 15 ml, ID</b>																				
Item	30314832																			
For	FC5718	FC5718R			▲		▲	▲		●										
Max Speed	4500 rpm	4500 rpm																		
Max RCF	2830 × g	2830 × g																		

\* Blood collection tube  
ID: auto rotor identification  
Hi: version with higher speed

● Adapter not required

▲ With adapter  
(Continue on next page)

<b>FC5718/FC5718R</b>		0.2/ 0.4 ml	0.5 ml	1.5/ 2.0 ml	4.5 ml	5 ml	7 ml	10 ml	15 ml	16 ml	25 ml	30 ml	50 ml	85 ml	100 ml	4- 10 ml *	5- 7 ml *	PCR	MTP
<b>Angle rotor, 20 x 10 ml, Hi</b>																			
Item	30314833																		
For	FC5718 FC5718R							●											
Max Speed	12000 rpm																		
Max RCF	15775 x g 15775 x g																		
<b>Angle rotor, 12 x 15 ml, ID</b>																			
Item	30314834																		
For	FC5718 FC5718R			▲			▲		●										
Max Speed	6000 rpm																		
Max RCF	4427 x g 4427 x g																		
<b>Angle rotor, 44 x 1.5/2.0 ml, ID</b>																			
Item	30314835 (30372718)																		
For	FC5718 FC5718R		▲						●										
Max Speed	13500 rpm																		
Max RCF	17113 x g 17113 x g																		
<b>Angle rotor, 30 x 1.5/2.0 ml, ID, Sealable</b>																			
Item	30314836																		
For	FC5718 FC5718R		▲						●										
Max Speed	13000 rpm																		
Max RCF	17758 x g 20595 x g																		
<b>Angle rotor, 24 x 1.5/2.0 ml, ID</b>																			
Item	30314838																		
For	FC5718 FC5718R		▲						●										
Max Speed	14000 rpm																		
Max RCF	18624 x g 21379 x g																		
<b>Angle Rotor, 24x1.5/2.0 ml, ID, BIOSEALS</b>																			
Item	30304361																		
For	FC5718 FC5718R		▲						●										
Max Speed	14000 rpm																		
Max RCF	18624 x g 21379 x g																		
<b>Angle rotor, 12 x 1.5/2.0 ml, ID</b>																			
Item	30314839																		
For	FC5718 FC5718R		▲						●										
Max Speed	18000 rpm																		
Max RCF	23542 x g 23542 x g																		
<b>Angle rotor, 64 x 0.5 ml, ID</b>																			
Item	30314840																		
For	FC5718 FC5718R																		
Max Speed	13500 rpm																		
Max RCF	16298 x g 16298 x g																		
<b>Angle rotor, 4 x 8-w PCR Strip, ID</b>																			
Item	30314841																		
For	FC5718 FC5718R																		
Max Speed	15000 rpm																		
Max RCF	15343 x g 15343 x g																		

Distributed by:  
 Safety Emporium  
 PO Box 1003  
 Blackwood, NJ 08012  
 Ph: (866) 326-5412 toll-free  
 Fax: (856) 553-6154  
 esupport@safetyemporium.com  
 www.safetyemporium.com

\* Blood collection tube  
 ID: auto rotor identification  
 Hi: version with higher speed

● Adapter not required      ▲ With adapter